

STANDOUT

BLOOD AND PROTEIN REMOVER

PRODUCT USE:

STANDOUT is a quick remover of blood, protein, albumin etc., with maximum safety to fabrics. A buffered product, **STANDOUT** removes difficult stains. Easy to use; merely squirt on, work in, flush out. Penetrates rapidly and rinses out.

ADVANTAGES:

- Low use cost.
- Removes broad range of stains.
- Buffered formula prevents setting of stains.
- Rapid, sure stain removal with maximum safety.
- Stable formula for longer shelf life.

PRODUCT DESCRIPTION:

Form: Colourless, clear viscous liquid
pH: 10.5 - 11.5
S.G.: 0.99 ±.01

DIRECTIONS FOR USE:

Lightly scrape or dry-brush stains to remove all dried garments. Moisten stain thoroughly with water or steam gun, and apply a small amount of **STANDOUT**. Work the stain with a spatula until the stain is dissolved. Flush with clean water or steam.

STANDOUT should be used with caution on wool and silk garments.

On questionable fabrics, **STANDOUT** can be diluted with Mr. Signal (1 part **STANDOUT**, 2 parts Mr. Signal).

Ref. No. 11/17

PRECAUTIONS:

Keep out of reach of children.

FIRST AID:

Eyes: Flush with plenty of water.
Internal: Drink large amounts of water. Get medical attention!

Every endeavor has been made to ensure that the information in this Data Sheet is reliable, but we do not accept any liability for any loss, injury or damage which may result from its use.